AT53077

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Thank you for purchasing our product, please read the instructions for use carefully. Do not use the product before reading and understanding this manual. If you do not understand the warnings, cautions and recommendations, please contact your health care professional or dealer to avoid product damage and personal injury.

CAUTION: Inspect all parts of the product for damage caused in transit. If you notice such damage, DO NOT use the product. Contact the manufacturer for more information.

DESCRIPTION AND FEATURES

The corrective toe brace is a specially designed orthopedic device designed to correct and stabilize a toe that has deformed laterally (toward the other toes). The device is designed to be soft and comfortable even when worn with shoes. It is made of high-quality elastic neoprene, and provides gentle straightening of the big toe, which is in a valgus position.

INDICATIONS: Physical or mental limitations (e.g., visual impairment) that prevent safe handling of the product.

- Severe foot deformities: If the toe deformity is very advanced, it may be difficult or impossible to achieve effective correction with the device.
- Inflammatory conditions or infections: If there are inflammatory conditions, such as arthritis, or skin infections on the foot, wearing braces can worsen symptoms and delay healing.
- Skin prone to ulceration: If the skin on the foot is very sensitive or prone to ulceration, wearing braces can cause chafing or excessive pressure, which can lead to more serious problems.
- Blood circulation disorders: People with blood circulation disorders, such as diabetes or atherosclerosis, may be limited in their ability to wear braces due to the risk of compromised circulation and the possibility of ulcer formation.
- Foot abnormalities requiring surgical intervention: In cases where the toe deformity is very advanced and requires surgical intervention, wearing a brace may be ineffective or even harmful.
- Hallux valgus of the toe (hallux valgus)
- Pain and discomfort associated with hallux valgus of the toe
- Prevention of further progression of the deformity
- Improvement of foot function
- Rehabilitation after surgery
- Prevention of additional injuries

APPLICATION

The use of a corrective toe brace is intended to halt the progressive deformity, reduce pain associated with toe valgus, and restore the natural position of the joint. The orthosis can be worn inside footwear, allowing the patient to maintain an active lifestyle while supporting the healing process.

SIZE: universal

NOTES:

In the event of a device-related „serious incident“ that directly or indirectly led to, may have led to, or may lead to any of the following:

- (a) death of a patient, user or other person, or
- (b) temporary or permanent deterioration of the health of a patient, user or other person, or
- (c) a serious threat to public health

the above „serious incident“ must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient resides. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

NOTES:

In case of pain, allergic reactions or other disturbing symptoms unclear to the user related to the use of the medical device, consult a health care professional.

NOTE: It is forbidden to use the product other than in accordance with its intended use.

NOTE: The manufacturer is not responsible for damage caused by negligent cleaning and maintenance, or resulting from failure to follow the recommendations in this manual.

TARGET PATIENT GROUP


Persons struggling with diseases, dysfunctions or injuries, for the treatment/rehabilitation or compensation of which this device is intended (see the section on the purpose of the device of this manual). The device may be purchased by the user independently or on the recommendation of a physician, therapist or other specialist. Whether purchasing the device on its own or as indicated by a doctor/therapist/other specialist, it is important to consider the available sizes/necessary functions/sizes/variants of the device, indications and contraindications for use, as well as information provided by the manufacturer.

USE


Follow the steps below to put on the product.

Note: the product should be worn on a dry, clean, bare foot.


1. Unfasten the Velcro, located by the tags.

2. Insert the toe into the loop, so that the visible point of intersection of the strap is at the top of the toe.

3. Wrap the longer strap over the ankle and connect to the shorter strap by fastening the Velcro in such a way that the corrective toe brace does not put too much pressure on the leg or foot.

4. Check the comfort of use by walking. If there is pain or discomfort, try adjusting the Velcro fastening individually.

CLEANING AND MAINTENANCE

The product is recommended to be washed by hand with mild detergent at up to 40°C. Do not wash in a washing machine. Do not use strong detergents, bleach or rough objects for cleaning. Dry flat laid out in the shade, do not expose to direct sunlight. Do not tumble dry. Do not dry-clean. Do not iron.

	Hand wash at up to 40°C
	Do not bleach
	Do not iron.
	Do not tumble dry
	Do not dry clean

STORAGE

Store the product in a dry cool room, away from moisture and direct sunlight.

HOW TO DISPOSE OF THE PRODUCT

After the device is taken out of service, the medical device can be disposed of as ordinary municipal waste

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WARRANTY CARD

Model:
Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawisłańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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