AT53048
INSTRUKCJA UŻYTKOWANIA ORTEZA KRĘGOSŁUPA
USER MANUAL FOR SPINAL ORTHOSIS
BENUTZERHANDBUCH FÜR WIRBELSÄULENORTHESSEN
UŽIVATELSKÁ PŘÍRUČKA PRO ORTÉZU ZAD
NÁVOD NA POUŽITIE – ORTÉZA CHRBTICE
MANUEL D'UTILISATION DES ORTHÈSES RACHIDIENNES
GEbruikersHANDLEIDING VOOR SPINALE ORTHESEN
MANUAL DE INSTRUCCIONES PARA ORTESIS DE COLUMNA
MANUALE D'USO DELLE ORTESI SPINALI
ANVÄNDARMANUAL FÖR SPINALORTOSER

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EN

Thank you for purchasing our product, please read the instructions for use carefully. Do not use the product before reading and understanding this manual. If you do not understand the warnings, cautions and recommendations, please contact your health care professional or dealer to avoid product damage or personal injury.

CAUTION: Inspect all parts for damage caused in transit. If you notice such damage, DO NOT use the product. Contact the manufacturer for more information.

APPLICATION/INDICATIONS

The Spinal Orthosis with Pulley System is crafted from high-quality materials, ensuring excellent durability, comfort, and a proper fit. This orthosis is ideal for recovery after spinal surgery and for providing spinal support during recreational sports. A mechanical pulley system facilitates effective abdominal compression in various situations, delivering strong and consistent spinal support. Designed for both daily and sports use, the product ensures maximum support while maintaining comfort. Immediate medical consultation is advised if the orthosis is used post-injury. If the injury is severe or if there are any doubts, seek medical advice without delay. Adult supervision is required for children's use. Keep this product out of reach of small children.

- lower back pain syndrome
- degenerative disease of the spine
- discopathy of the lumbar segment

CONTRAINDICATIONS

Physical or mental limitations (e.g., visual impairment) that prevent safe handling.

- Use of the orthosis on an area with open wounds or skin infections may aggravate wounds and increase the risk of spreading infection.
- An orthosis that does not fit properly can lead to additional problems, such as abrasions, pain or worsening existing spinal problems.
- In some cases, advanced neurological diseases may preclude the use of a spinal orthosis, as the additional stiffness may adversely affect the patient's mobility or worsen other neurological symptoms.
- If the patient is unable to put on and take off the orthosis on his or her own, and no assistance is available, this may be a contraindication to its use.
- Patients with sensory disturbances in the area covered by the orthosis run the risk of not feeling abnormal pressure or other problems that can lead to skin or soft tissue injuries.
- In some cases, especially with orthoses involving the upper spine (such as thoracolumbar orthoses), the use of the orthosis may limit the ability to breathe deeply, which is particularly dangerous for those with pulmonary problems.

FEATURES

- Stiff and stable 22-centimeter lumbar panels direct stabilizing force along the lumbar vertebrae(L1-5), protecting against excessive lateral and rotational movement
- Double tightening system allows separate adjustment - the right tightener is responsible for the right part of the panels, and the left tightener for the left part of the panels
- the belt, made of free-breathing mesh, has a lightweight, open structure, promoting the removal of heat and sweat through holes in the material
- thanks to the pulley system, the orthosis has a wide range of adjustments
- free control of support and comfort; an easy-to-tighten drawstring with an opening in the middle, allows for optimal compression

SIZE

TO FIT THE RIGHT SIZE OF THE ORTHOSIS, MEASURE THE WAIST CIRCUMFERENCE.

S	M	L	XL	XXL
56-71 CM	69-84 CM	81-96 CM	94-109 CM	107-122 CM

HEIGHT S,M,L,X-22CM,XXL-22.5CM

COMPOSITION: 44% polyethylene, 31% nylon, 25% dacron

COMPONENTS:

Spinal orthosis, Instructions for use.

NOTES:

In the event of a device-related "serious incident" that directly or indirectly led to, may have led to, or may lead to any of the following:

- (a) death of a patient, user or other person, or
- (b) temporary or permanent deterioration of the health of a patient, user or other person, or
- (c) a serious threat to public health

the above "serious incident" must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient resides. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

CAUTION:

Before each use, make sure that the orthosis is free of seam abrasions, material tears, etc. Use of the product contrary to the rules contained in the instructions may result in serious complications. The orthosis should be used as recommended by a doctor or physiotherapist. -The duration of wearing the orthosis should be in accordance with the recommendations of the doctor or physiotherapist. The product should not be used in case of open wounds, abrasions of the epidermis, or after previous application of warming gel or ointment.

CAUTION:

In case of pain, allergic reactions or other disturbing, unclear to the user symptoms related to the use of the medical device, consult a healthcare professional.

NOTE: It is forbidden to use the product other than in accordance with its intended use.

NOTE: The manufacturer is not responsible for damage caused by negligent cleaning and maintenance, or resulting from failure to follow the recommendations in this manual.

TARGET PATIENT GROUP

Persons struggling with diseases, dysfunctions or injuries, for the treatment/rehabilitation or compensation of which this product is intended (see the section on the purpose of the product of this user manual). The device may be purchased by the user independently or on the recommendation of a physician, therapist or other specialist. Whether purchasing the device on its own or as indicated by a doctor/therapist/other specialist, the available sizes/necessary functions and variants of the device, indications and contraindications for use, as well as information provided by the manufacturer should be taken into account.

USE

1. The orthosis should be worn in a standing position.
2. Place the pulley system in the middle of the back.
3. Place the left side of the orthosis on your abdomen, then attach the right side to it with Velcro. Make sure the orthosis is compressed enough so that it does not fall off.
4. Unfasten the side hooks with holes in the middle, then pull forward which will compress the orthosis properly.

5. After you get the proper support(no pain, pressure, discomfort, feel comfortable and free) attach the handles with holes to the Velcro.

NOTE: Unless otherwise instructed by a doctor, the orthosis should be used intermittently, for a total of a maximum of 6 hours per day. Failure to do so may result in relaxation of the peroneal muscles. Unless otherwise advised by your doctor, do not use the orthosis while sleeping.

CLEANING AND CARE

The orthosis can only be washed by hand in lukewarm water (up to 30°C) with a mild detergent (such as gray soap).

- Do not iron. - Do not dry wash. - Do not bleach. - Protect from moisture. - Do not expose to high temperatures. - Do not tumble dry. - Dry unfolded, squeeze gently before drying.

STORAGE

Store the product in a dry cool room, away from moisture and direct sunlight.

HOW TO DISPOSE OF THE PRODUCT

After the product is taken out of service, the medical device can be disposed of as ordinary municipal waste.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiślańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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