AT53111

**INSTRUKCJA UŻYTKOWNIA, ORTEZA KRĘGOSŁUPA
USER MANUAL, SPINAL ORTHOSIS
GEBRAUCHSANWEISUNG, WIRBELSÄULENORTHESSE
NÁVOD K POUŽITÍ, ORTÉZA PÁTEŘE
NÁVOD NA POUŽITIE, ORTÉZA CHRBTICE
MODE D'EMPLOI, ORTHÈSE SPINALE
GEBRUIKSAANWIJZING, WERVELKOLOMORTHESE
INSTRUCCIONES DE USO, ÓRTESIS DE COLUMNA
ISTRUZIONI PER L'USO, ORTESI SPINALE
BRUKSANVISNING, RYGGRADSORTOS**

MD

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En

Thank you for purchasing our product. Please read this manual carefully before use. Do not use the product before you have read and understood this instruction manual. If you do not understand the warnings, notes, or recommendations, please contact a healthcare professional or the distributor to avoid product damage or personal injury. **WARNING:** Check all parts for any damage that may have occurred during transport. If such damage is found, DO NOT use the product. For more information, contact the seller.

INTENDED USER GROUP

Individuals with diseases, dysfunctions, or injuries for whose treatment, rehabilitation, or compensation this product is intended (see the section "Purpose of the Device" in this manual).

The product can be purchased independently by the user or on the recommendation of a doctor, therapist, or other healthcare specialist. When purchasing independently or based on medical recommendation, consider the available sizes, required functions, product variants, indications and contraindications for use, as well as the information provided by the manufacturer.

APPLICATION

The spinal brace is a medical device designed to stabilize, support, or correct the spine.

It is intended for the treatment of injuries, conditions, or dysfunctions of the musculoskeletal system in the spinal region.

The brace limits spinal mobility, provides support, and relieves the skeletal, muscular, and ligamentous structures after injuries, surgeries, or in cases of strain.

INDICATIONS

Injuries and fractures of the lumbar spine requiring external stabilization.

Overload, instability, or weakening of muscles in the lumbar-sacral region.

Lumbar disc disease and other conditions causing lower back pain.

Post-traumatic and post-operative states in the lumbar spine area.

Preventive use to protect the spine in individuals performing physical work or spending long periods in standing or sitting positions.

CONTRAINDICATIONS

Physical or mental limitations (e.g., visual impairment) preventing safe handling of the product. Acute skin inflammation, open wounds, or ulcers in areas in contact with the brace. Known allergy or hypersensitivity to any materials used in the product (e.g., latex, neoprene).

Skin lesions, rashes, or irritation in the lumbar region that may worsen due to pressure. Severe spinal deformities preventing proper fit of the brace. Use without medical consultation in cases of recent fractures, injuries, or spinal tumors.

WARNING: The brace contains strong magnets. Do not use in persons with a pacemaker, ICD, or other implanted electrical devices. The magnetic field may interfere with their function and pose a life-threatening risk.

COMPONENTS: Spinal brace, User manual

TECHNICAL SPECIFICATION

Sizes:

M – Length: 90 cm, Suitable waist circumference: 65–80 cm

L – Length: 100 cm, Suitable waist circumference: 75–95 cm

XL – Length: 110 cm, Suitable waist circumference: 85–105 cm

2XL – Length: 120 cm, Suitable waist circumference: 95–115 cm

3XL – Length: 130 cm, Suitable waist circumference: 105–125 cm

cm

Back height: 24.5 cm Front height: 16 cm

USE INSTRUCTIONS

Material composition:

Main fabric: 40% polyester / 10% nylon / 35% neoprene / 15% latex

Metal splint: steel and aluminum alloy plate, PP plate

Velcro: nylon

Strap: polyester, latex

Wrap the brace around the lower back so that the support bars are centered along the spine.

Pull the main straps forward and fasten them with Velcro.

Tighten the side straps for proper fit and support.

The brace should fit snugly but not cause discomfort.

If pain or discomfort occurs, loosen the straps or remove the product.

WARNINGS

In the event of a serious incident related to the device that has led, might have led, or could lead to any of the following:

a) death of a patient, user, or other person,

b) temporary or permanent deterioration in the health of a patient, user, or other person, or

c) a serious threat to public health,

such an incident must be reported to the manufacturer and the competent authority of the Member State where the user or patient resides. **NOTE:** If pain, allergic reactions, or any unusual or concerning symptoms occur during use, discontinue use and consult a healthcare professional.

WARNING: The product must not be used for purposes other than those for which it is intended.

CAUTION: This product contains natural latex, which may cause allergic reactions in individuals sensitive to latex.

CLEANING AND MAINTENANCE: Hand wash in cold water with mild soap. Air-dry away from heat sources and sunlight. Do not bleach, iron, dry-clean, or tumble dry. **STORAGE:** Store in a dry place, away from direct sunlight. **DISPOSAL** After the end of its service life, dispose of the product as standard household waste. Follow local waste segregation regulations — protect the environment.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR MEDICAL Sp. z o.o., 03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR MEDICAL Sp. z o.o.
03-068 Warsaw, ul. Zawisłańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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