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INSTRUCTION MANUAL
RIGID TORSO ORTHOSIS AT04501
LUMBOSACRAL ORTHOSIS AT04502, AT04506
LUMBOSACRAL OR PELVIC STABILIZING BELT AT04503
RIGID TORSO ORTHOSIS HIGH WITH STACKS AT04505
STABILIZATION AND IMMOBILIZATION CORSET AT04510

Thank you for purchasing our product, please read the instructions for use carefully. Do not use the product before reading and understanding this manual. If you do not understand the warnings, cautions and recommendations, please contact your health care professional or dealer to avoid damage to the product.

NOTE: Inspect all parts for damage caused in transit. If you notice such damage, DO NOT use the product. Contact the manufacturer for more information.

APPLICATION AT04502/3/6/10

The spinal orthosis improves stabilization of the lumbar spine, reducing pain and improving the patient's daily comfort. Stabilization is provided by elastic straps and contoured reinforcements (nibs), which relieve pressure on the spinal tissues by transferring forces that weigh down the spinal structures, preventing harmful changes in the lower spine.

APPLICATION AT04501/5

The spinal orthosis improves stabilization of the lumbar and thoracic spine reducing pain and improving the patient's daily comfort. Stabilization is ensured by elastic fastening straps, additional reinforcements and straight nibs, which relieve the pressure on the spinal tissues by transferring the forces weighing down the spinal structures.

INDICATIONS

- spinal disorders
- multilevel discopathy
- osteoporosis
- degenerative changes
- spinal instability
- lumbar chiropractic,
- Osteoporotic compression fracture of the spine
- sciatica
- lumbosacral pain syndromes
- period after surgical procedures

CONTRAINDICATIONS

physical or mental limitations (e.g., visual impairment) that prevent safe handling.

CAUTION: Do not use the product after prior use of warming gels and ointments.

FEATURES AT04502

- The orthosis has 6 steel reinforcements: 4 with a length of 28cm and 2 with a length of 24cm)

FEATURES AT04503

- The orthosis has 4 steel reinforcements with a length of 20 cm

FEATURES AT04506

- The orthosis has 4 steel reinforcements with a length of 20 cm

FEATURES AT04510

- The orthosis has 4 steel reinforcements with a length of 28 cm

FEATURES AT04501

- The orthosis has 2 steel reinforcements with a length of 44 cm

FEATURES AT04505

- The orthosis has 4 steel reinforcements with a length of 28 cm and 2 straight reinforcements with a length of 44 cm

SIZES AT04502

S (80 - 90 cm), M (90 - 100 cm), L (100 - 110 cm), XL (110 - 120 cm),
XXL (120 - 130 cm), XXXL (130 - 140 cm), XXXXL (140 - 150 cm)

SIZES AT04503

XS (70 - 80 cm), S (80 - 95 cm), M (90 - 100 cm), L (95 - 105 cm),
XL (105 - 115 cm), XXL (115 - 125 cm), XXXL (125 - 135 cm)

SIZES AT04506

S (60 - 75cm), M (76 - 90 cm), L (91 - 100 cm), XL (101 - 115 cm),
XXL (116 - 125 cm), XXXL (126 - 140 cm)

SIZES AT04510

S (80-95cm), M (90-100cm), L (95-105cm), XL (105-115cm),
XXL (115 - 125 cm), XXXL (125 - 135 cm)

SIZE AT04501 (To choose the right size, measure your waist circumference)

S (85-95 cm), M (95-105 cm), L (105-115 cm), XL (115-125 cm),
XXL (125-136 cm), XXXL (136-145 cm)

SIZES AT04505

S (80 - 95cm), M (90-100cm), L (95-105cm), XL (105-115cm), XXL (115 -125
cm), XXXL (125-135cm)

HEIGHT:

AT04502 - 30 cm

AT04503 - 23 cm

AT04506- 24 cm

AT04510 - 30 cm

AT04501- 54 cm

AT04505 - 54 cm

COMPONENTS

Spinal orthosis, instruction manual

NOTES:

In the event of a product-related "serious incident" that directly or indirectly led to, may have led to, or may lead to any of the following:

- (a) death of a patient, user or other person or
- (b) temporary or permanent deterioration of the health of the patient, user or other person or
- (c) a serious threat to public health

the above "serious incident" should be reported to the manufacturer and to the competent authority of the Member State in which the user or patient resides. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

NOTES:

If you experience pain, allergic reactions, or other disturbing symptoms that are unclear to you related to the use of the medical device, consult a healthcare professional.

NOTE: it is forbidden to use the product other than for its intended purpose .

NOTE: The manufacturer is not responsible for damage caused by negligent cleaning and maintenance, or resulting from failure to follow the recommendations in this manual.

TARGET PATIENT GROUP

Persons struggling with diseases, dysfunctions or injuries, for the treatment/rehabilitation or compensation of which this device is intended (see the section on the purpose of the device of this manual). The device may be purchased by the user independently or on the recommendation of a physician, therapist or other specialist. Whether purchasing the device on its own or as indicated by a doctor/therapist/other specialist, the available sizes/necessary functions and variants of the device, indications and contraindications for use, as well as information provided by the manufacturer, should be taken into account.


USE AT04502/3/6/10

1 Make sure the orthosis is in the correct position. Put on the orthosis by fastening the main strap with Velcro, tighten the abdominal strap tightly enough.

(2) Then fasten the two additional tightening straps by grasping both ends at the same time.


3. Finally, once again make sure that the device properly adheres to the spine, is put on at the right height, the straps are squeezed tightly enough and provide stabilization.

NOTE: We recommend that a doctor or specialist be present when trying to use the product for the first time.

1 Make sure that the lower belt of the product is at the level of the lumbar spine and the additional long belts are on the shoulders. Put on the orthosis by fastening the main belt with Velcro, be sure to squeeze the abdominal belt tightly enough. (AT045010) 2. Then fasten the two additional tightening straps by grasping both ends at the same time. (AT045010) 3. pass the additional long straps under the shoulders and cross behind the back, pull forward and fasten on the main belt. 4. Finally, once again make sure that the device properly adheres to the spine, is put on at the right height, the straps are squeezed tightly enough and provide stabilization.

NOTE: We recommend that a doctor or specialist be present when trying to use the product for the first time.

CLEANING

The product is recommended to be washed by hand with a mild detergent, at a temperature of up to 40 ° C. Do not wash in a washing machine. Do not use strong detergents, bleach or rough objects for cleaning. Dry flat laid out in the shade, do not expose to direct sunlight. Do not tumble dry. Do not dry-clean. Do not iron.

STORAGE

Store the product in a dry cool room, away from moisture and direct sunlight.

METHOD OF DISPOSAL OF THE PRODUCT

The medical device can be disposed of like ordinary municipal waste.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiślańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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