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AT04202

AT53078

AT04203 Ankle stabilizing orthosis (with cold warm insert)


AT04204 Ankle stabilizing orthosis (with foam insert and pneumatic insert)

AT53078 Ankle stabilizing orthosis (with cold warm insert)

The orthosis provides excellent stiffening of the ankle joint, without excessively limiting the mobility of the foot. The stabilizers are recommended during rehabilitation after 2nd and 3rd degree sprains of the ankle and foot joint,

Use of AT04201

1. Fix the inserts inside the plastic covers.
2. Adjust the position of the heel pad with the side straps.
3. Position the stabilizer so that the cube is in the recesses intended for it.
4. Attach the tightening straps.

Insole length: 25 cm

Insole width: 10 cm at widest point Velcro length: 46 cm

Recommendations:

Before applying the stabilizer, wear a thick sock;

The stabilizer should be cleaned with water and possibly soap. The insoles should be washed by hand.

Do not use detergents. The insoles should not be ironed.

AT53078

Ankle orthosis with internal gel pad improves circulation and reduces swelling. It can be heated or cooled during thermotherapy. Stiff, anatomically designed guards provide protection, comfort and prevent ankle sprain. Velcro closures provide a comfortable and easy fit.


INDICATIONS/APPLICATION

- Sprains of the ankle joint
- Tendonitis
- Ligament injuries
- Chronic instability of the ankle joint
- Ankle sensitivity

SIZE: UNIVERSAL USE


The orthosis should be put on according to the following order.

1. Place the orthosis on the ankle joint.
2. Adjust the position of the heel pad with the side straps.
3. Position the stabilizer so that the ankle is in the recesses intended for it.
4. Attach the tightening straps.

AT04202 stabilizer application:

1. Position the stabilizer so that the ankle is in the recesses intended for it.
2. Fasten the tightening straps.
3. Place the end of the plastic hose in the valve located in the top of the insert.
4. Orally inflate the orthotic insert. Pump until the ankle properly stabilized.
5. After inflating, remove the hose and tighten the valve with your fingers. Proceed in the same way with the second insert.
6. To drain the air from the insert, put the hose back in the valve and gently press on the insert to drain the air.

TARGET PATIENT GROUP

Persons struggling with diseases, dysfunctions or injuries, for the treatment/rehabilitation or compensation of which this product is intended (see the section on the purpose of the product of this manual).

The device can be purchased by the user independently or based on the recommendation of a doctor, therapist or other specialist. Whether purchasing the device on its own or as indicated by a doctor/therapist/other specialist, it is important to consider the available sizes/necessary functions and variants of the device, indications and contraindications for use, as well as information provided by the manufacturer.

CLEANING

Clean the product with a damp cloth with mild detergent and then wipe with a dry cloth. Do not use rough objects, solvents or strong detergents.

STORAGE

Store the product in a dry cool room, away from moisture and direct sunlight.

DISPOSAL OF THE PRODUCT

Dispose of the product in accordance with the applicable disposal conditions.

CAUTION: it is forbidden to use the product other than for its intended purpose.

NOTE: The manufacturer is not responsible for damage caused by negligent maintenance, inadequate servicing or resulting from failure to follow recommendations in this manual.

NOTES:

In the event of a product-related "serious " that directly or indirectly led to, may have led to, or may lead to any of the following:

- a) death of a patient, user or other person or
- b) temporary or permanent deterioration in the health of a patient, user or other person or
- c) a serious threat to public health

the above "serious incident" should be reported to the manufacturer and to the competent authority of the Member State in which the user or patient resides. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

NOTES:

If you experience pain, allergic reactions, or other disturbing symptoms that are unclear to you related to the use of the medical device, you should stop using the device and consult a healthcare professional.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiślańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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