AT51129

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USER MANUAL

AT51129 BACK SUPPORT

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AT51129 OPĚRKA ZAD

NÁVOD NA POUŽITIE

AT51129 OPIERKA CHRBT

MODE D'EMPLOI

AT51129 SOUTIEN DU DOS

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AT51129 RUGSTEUN

MANUAL DE INSTRUCCIONES

AT51129 APOYO LUMBAR

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AT51129 SUPPORTO SCHIENALE

BRUKSANVISNING

AT51129 RYGGSTÖD

PL

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EN

Thank you for purchasing our product. Please carefully read the user manual. Do not use the product before reading and understanding this manual. If you do not understand the warnings, notes, and recommendations, contact a healthcare specialist or the seller to avoid damaging the product or causing personal injury.

WARNING: Check all parts for damage caused during transport. If such damage is noticed, DO NOT use the product. More information is available from the seller.

TARGET PATIENT GROUP

People suffering from diseases, dysfunctions, or injuries for which this product is intended for treatment/rehabilitation or compensation (see the section on the intended use of this product in this user manual). The product may be purchased by the user independently or based on recommendations from a doctor, therapist, or other specialist. Whether purchased independently or recommended by a doctor/therapist/specialist, the available sizes/necessary functions and product variants, indications, and contraindications for use, as well as the information provided by the manufacturer, should be considered.

USE/INDICATIONS

The support is designed to allow bedridden patients to sit up in bed for taking medications and eating, as well as to transition from a lying to a sitting position. Such a change in position helps improve blood circulation in bedridden individuals, increasing their comfort and supporting regenerative processes.

People with limited mobility: Patients who have difficulty maintaining a sitting position due to muscle weakness, paralysis, neurological diseases, or other conditions may require back support.

Respiratory diseases: Patients with respiratory problems, such as COPD, asthma, or heart failure, benefit from sitting with appropriate back support, which helps improve breathing by reducing pressure on the lungs and facilitating ventilation.

Post-operative patients: Patients who have undergone surgical procedures, particularly in the abdominal or chest area, may need support to avoid pain and discomfort while sitting and eating.

Pressure ulcer prevention: For chronically bedridden patients, proper back support can prevent the development of pressure ulcers by ensuring even pressure distribution on the back.

Neurological diseases: For individuals with neurological diseases, such as stroke, multiple sclerosis, or Parkinson's disease, back support can help maintain a stable and safe sitting position during meals.

Dysphagia (swallowing difficulties)

Elderly individuals

CONTRAINDICATIONS

Physical or mental limitations (e.g., visual impairment) that prevent safe handling of the product.

For patients with spinal injuries, particularly in the cervical or thoracic regions, raising to a sitting position may be dangerous. In such cases, specialist recommendations or orthopedic cushions should be used.

When a patient experiences severe back, neck, or hip pain, changing positions may exacerbate symptoms. In such cases, consultation with a doctor or physiotherapist is necessary to adjust the position in a way that minimizes discomfort.

Patients with serious cardiac issues, such as orthostatic hypotension (sudden drops in blood pressure upon changing position), may experience dizziness, weakness, or fainting when raised to a sitting position.

Patients with chest injuries, such as fractured ribs or pneumothorax, may feel discomfort or pain when attempting to sit up. Special adjustments to the support are necessary in such situations.

In some cases, patients with advanced respiratory failure may not tolerate the sitting position if it increases pressure on the diaphragm, worsening respiratory function. In these cases, the height of the support should be adjusted, or the full sitting position should be avoided.

In some situations, especially post-abdominal surgeries, too high a sitting position may cause tension in the operated area, delaying healing or causing pain.

Patients with neurological disorders that affect movement control or balance (e.g., severe stroke, amyotrophic lateral sclerosis) may be at increased risk of falling or tipping sideways in an unsupported sitting position.

FEATURES

Easy to use

Backrest adjustment (6 levels)

Headrest pillow

Breathable fabric (100% polyester)

Steel Frame

TECHNICAL SPECIFICATIONS

Folded size: 60x54.5x8cm

Maximum size when unfolded: 56x54.5x50cm

Weight: 1.75kg

Maximum load: 150kg

This symbol indicates the maximum user weight.

USAGE INSTRUCTIONS

Unfold the support and place it on the bed in a position that will support the patient's back.

Adjust the backrest angle to the patient's preference using the available adjustments on the support.

Position the patient to ensure comfortable support for the back and head.

Note: Ensure the patient's hands are away from the adjustment mechanisms when changing the backrest angle.

Note: Verify that the locking mechanism for the angle adjustment is correctly engaged before use to prevent sudden collapse.

Note: Ensure the product is stable and securely positioned before each use.

COMPONENTS

Back support, user manual

WARNING:

In the event of a "serious incident" related to the product that directly or indirectly resulted, could have resulted, or may result in any of the following events:

- a) the death of the patient, user, or another person, or
- b) temporary or permanent deterioration of the health of the patient, user, or another person, or
- c) a serious public health threat,

the "serious incident" must be reported to the manufacturer and the competent authority of the member state where the user or patient resides. In Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices, and Biocidal Products.

WARNING:

If pain, allergic reactions, or other concerning or unclear symptoms related to the use of the medical device occur, consult a healthcare professional.

WARNING:

It is prohibited to use the product in a manner other than as intended.

WARNING: If the patient cannot move independently, their position should be changed regularly. Otherwise, pressure sores may develop.

CLEANING AND MAINTENANCE

Use a mild household cleaner to clean metal components and, if necessary, a standard disinfectant. Do not use aggressive, corrosive, or abrasive cleaners.

Regularly check for missing parts and ensure that all components are in good technical condition. If wear, mechanical damage, or missing parts are noticed, discontinue use immediately and contact the manufacturer's service department to avoid the risk of accidents or further damage.

Use a damp cloth to clean fabric components.

STORAGE

Store the product in a dry place and do not expose it to direct sunlight.

DISPOSAL

Once the medical device has been decommissioned, it can be disposed of as regular municipal waste. For proper disposal, please contact your local authority or nearest waste collection center to learn how to properly dispose of the product. By complying with local regulations, you help protect the environment.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiślańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

 ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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