AT51134
INSTRUKCJA UŻYTKOWNIA, CHODZIK
INSTRUCTION MANUAL, ROLLATOR
GEBRAUCHSANWEISUNG, ROLLATOR
NÁVOD K POUŽITÍ, CHODÍTKO
NÁVOD NA POUŽITIE, CHODÍTKO
MANUEL D'INSTRUCTION, DÉAMBULATEUR
GEbruikSAANWIJZING, ROLLATOR
MANUAL DE INSTRUCCIONES, ANDADOR
MANUALE DI ISTRUZIONI, DEAMBULATORE
BRUKSANVISNING, RULLATOR

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EN

Thank you for purchasing our product. Please read this user manual carefully. Do not use the product before you have familiarized yourself with and fully understood this manual. If you do not understand the warnings, notices, or recommendations, please contact a healthcare professional or the seller to avoid damage to the product or personal injury. **WARNING:** All parts must be checked for damage that may have occurred during transport. If any such damage is found, the product **MUST NOT** be used. For further information, please contact the seller.

TARGET PATIENT GROUP

People suffering from diseases, dysfunctions, or injuries for which this product is intended for treatment, rehabilitation, or compensation (see the section on intended use in this user manual). The product may be purchased independently by the user or based on the recommendation of a physician, therapist, or another specialist. Regardless of whether the product is purchased independently or recommended by a physician/therapist/other specialist, the available sizes, required functions, and product variants, as well as indications and contraindications for use and the information provided by the manufacturer, must be taken into account.

INTENDED USE / INDICATIONS

The product is intended to support mobility in persons with limited walking ability. The rehabilitation rollator is designed for use by one person only. The product may be used in the following cases:

Balance training: The rehabilitation rollator allows patients to practice maintaining balance, which is important in the rehabilitation process after injuries or neurological disorders.

Posture support: Patients can use the rehabilitation rollator to train correct body posture, which is particularly important for people with postural defects or problems maintaining body stability.

Improvement of motor coordination: Exercises performed using the rehabilitation rollator can help improve motor coordination and movement control, which is essential for people with neurological disorders.

Support during post-injury rehabilitation: The rehabilitation rollator can be used in the rehabilitation process after lower limb injuries, helping patients gradually regain mobility.

Increasing independence: For some patients, the rehabilitation rollator may be a tool to increase independence in everyday life by enabling independent movement.

CONTRAINDICATIONS

Physical or mental limitations (e.g. visual impairment) that prevent safe use of the product.

Individuals with serious spinal stability problems or spinal injuries may not be suitable candidates for using a rehabilitation rollator, as this may increase the risk of further injury.

In people with uncontrolled seizures or epilepsy, the use of a rehabilitation rollator may be dangerous both for the patient and for people nearby.

In individuals with very limited muscle strength, using a rehabilitation rollator may be difficult or ineffective. In such cases, alternative rehabilitation approaches may be necessary, such as manual therapy or muscle-strengthening exercises.

In patients with severe balance disorders, the use of a rollator may be unsafe, as it can increase the risk of falls and further injuries.

LIST OF COMPONENTS: Rollator, user manual

PRODUCT FEATURES

Ergonomic handles

Durable, foldable frame with locking mechanism

Adjustable bars

Brake system with a light and easy-to-use mechanism

Padded forearm support

Front and rear wheels suitable for difficult terrain and uneven surfaces

Damage-resistant construction suitable for home and institutional use

Foldable seat allowing easy access to a shopping bag or personal items

Frame designed to provide improved accessibility and freedom of movement while walking

Easy to transport and store

Maximum safe load: up to 136 kg

TECHNICAL SPECIFICATION

Frame: aluminium

Wheel fork: PA + 30% glass fiber

Chest support / support surface: foam and PVC, dimensions 60 × 50 × 6 cm

Handles: PVC 70%

Width between handles / distance between handles: 48 cm

Overall width: 64 cm

Overall depth: 86 cm

Chest support height (min–max): 108–131 cm

Adjustable handle height (min–max): 120–144 cm

Chest support thickness (support element): 6 cm

Folded dimensions (L × W × H): 69 × 26 × 110 cm

Front wheels: 8" (PVC + PP)

Rear wheels: 8" (PVC + PP)

Maximum user weight: 136 kg

Net weight: 11.3 kg

WARNING: Individual product dimensions may vary within a tolerance of ±2 cm.

This symbol indicates the maximum user weight.

ASSEMBLY

Unfold the rollator frame on a stable, flat surface.

Insert the black vertical tube into the opening located on both sides of the main frame.

Align the holes in the tube with the adjustment holes in the frame.

Make sure the tube is inserted to the correct depth.

Lock the height – press the locking pin into the adjustment hole until you hear a characteristic “click.”

Place the upper part with the handles onto the previously installed vertical tubes.

Tighten the fixing knobs located under the handles to securely attach the element to the frame. Tighten the knobs by hand until noticeable resistance is felt. Check that the entire assembly does not move when pressure is applied.

One side of the support must be secured using the screw and wrench included in the set. Make sure the screw is firmly tightened and the element is securely attached.

OPERATING INSTRUCTIONS

To brake while moving – squeeze the brake levers under the handles.

To lock the rollator in a stationary position – push the brake levers downward until they lock.

When sitting or standing still, the brakes must be engaged.

To fold the rollator – engage the brakes, press the black lever at the rear of the frame, and fold the base inward.

WARNINGS REGARDING THE USE OF THE ROLLATOR

Use of the rollator: The rollator must be used only as intended. Improper use may result in injury.

Stability: Always make sure that all legs/wheels of the rollator are in contact with the ground and stable before starting to walk. Walking on uneven or sloped surfaces may be dangerous.

Use on stairs: The rollator is not intended for use on stairs. Before going up or down stairs, use other aids or ask for assistance.

Maximum load: Do not exceed the maximum load specified by the manufacturer. Overloading the rollator may cause damage and increase the risk of accidents.

Transport: Always ensure that the rollator is properly folded or secured during transport.

Regular inspections: Before each use, check the rollator for damage, loose screws, or other irregularities. If any damage is found, stop using the rollator and contact service personnel.

Warning! The rollator is not a toy. Do not allow children to play with it or use it without adult supervision.

Be aware of your surroundings: Avoid using the rollator on slippery, wet, or uneven surfaces. Objects on the ground may cause tripping and falls.

IMPORTANT NOTES

When the rehabilitation rollator is not in use, make sure the parking brake is engaged.
Before use, check that all wheels are securely fastened and that the brakes function properly.
Check that the height adjustment button is locked.
Make sure the rollator is fully unfolded.
Pay attention to your surroundings while using the rollator.
Do not place your hands in the folding area to avoid pinching injuries.
Do not drop yourself onto the seat; always sit down slowly and carefully. Ask for assistance if necessary.
Do not hang any objects on the brake lever or handle. This may affect brake performance or rollator stability.
Do not stand on the rehabilitation rollator.
Make sure both handles are set to the same height.
WARNING: Use of the product in a manner other than its intended purpose is prohibited.

SERIOUS INCIDENTS

In the event of a "serious incident" related to the medical device that directly or indirectly has led, could have led, or may lead to any of the following events:

- a) death of the patient, user, or another person, or
 - b) temporary or permanent deterioration of the health of the patient, user, or another person, or
 - c) a serious public health threat,
- the serious incident must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient resides. **WARNING:** Use of the product in a manner other than its intended purpose is prohibited.

CLEANING

Clean the device regularly and keep it tidy.
Light dirt and dust can be removed with a damp cloth.
Heavier dirt can be cleaned with a mild, non-corrosive detergent and warm water.
Do not use strong abrasives or chlorine-based cleaners as they may cause damage or discoloration.
Dry parts with a soft cloth.
Wheels: Clean wheels with warm water and mild detergent. Dry with a clean cloth. Lubricate periodically with multi-purpose oil or when wheels start to squeak.

MAINTENANCE

To ensure safe and proper use, inspect the device regularly according to the schedule below:

Weekly checks:

- Check that all screws are tight.
- Inspect brake cables for looseness or damage.
- Tighten screws and nuts if necessary.

Monthly checks:

- Check handle stability and condition.
- Inspect brake levers for proper function.
- Inspect bases (legs) for stability and looseness.

Quarterly checks (every 4 months):

- Inspect wheels to ensure they rotate freely and are undamaged.
- Check the bag for wear or damage.
- Inspect the frame and supports for cracks or deformations.

TRANSPORT

The rollator folds easily, making it convenient to transport. For transportation, fold it as described above. Always transport the rollator securely to prevent sudden, unsafe movements.

REUSE

The rollator is suitable for reuse. Before passing it to another user, it must be inspected by a qualified dealer. Check the following:

- Visual condition of the frame.
- Secure fit of screws and nuts.
- Condition of plastic parts.
- Condition of wheels and bearings.
- Function and condition of brakes.

If any part is damaged or broken, it must be replaced. Wheels should be replaced whenever significant wear is noticed. All loosened screws and nuts must be tightened. The product must be thoroughly cleaned and disinfected before reuse.

WARNING: Do not reuse the rollator if the frame is bent or broken.

The user manual must always accompany the rollator when changing users.

STORAGE

Store the product in a dry place. Do not expose it to direct sunlight.

DISPOSAL

When taken out of service, dispose of the product as regular household waste. Follow local waste separation regulations. Protect the environment.

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Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR MEDICAL Sp. z o.o. ,03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.
10. The warranty does not cover parts and components subject to natural wear and tear, such as wheels, tires, rubber pads, fuses, etc., the wear of which results from normal use of the product. The warranty period for batteries is 6 months.

ANTAR MEDICAL Sp. z o.o.
03-068 Warsaw, ul. Zawisłańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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