AT51108
INSTRUKCJA UŻYTKOWNIA, LASKA DLA NIEWIDOMYCH
USER MANUAL, CANE FOR THE BLIND
GEBRAUCHSANWEISUNG, BLINDENSTOCK
NÁVOD K POUŽITÍ, HŮL PRO NEVIDOMÉ
NÁVOD NA POUŽITIE, PALICA PRE NEVIDIACICH
NOTICE D'UTILISATION, CANNE POUR AVEUGLES
GEbruiksaanwijzing, BLINDENSTOK
MANUAL DE INSTRUCCIONES, BASTÓN PARA PERSONAS CIEGAS
MANUALE DI ISTRUZIONI, BASTONE PER NON VEDENTI
BRUKSANVISNING, KÄPP FÖR BLINDA

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EN

READ THE ENTIRE INSTRUCTION MANUAL CAREFULLY BEFORE USE

PRODUCT DESCRIPTION:

Universal (can be used as both right- and left-hand version),
Made of aluminum,
Handle made of plastic,
Consists of four elements. Foldable,
Maximum load capacity: 120 kg,
Tip: black color,
Service life: 5 years from the date of manufacture,
Cane weight: 0.21 kg,
Length of the cane when unfolded: 120 cm,
Includes a carrying/storage pouch.

INTENDED USE:

The product is intended for a person with a disability, usually visually impaired. The cane helps the user independently orient themselves in space and warns of obstacles that may appear along the way. The cane performs orientation, warning, and protective functions. The cane consists of three aluminum sections connected with an elastic cord, finished with a tactile plastic tip, and features a handle with a wrist strap.

UNFOLDING:

The cane is supplied in a folded form. After removing it from the packaging, unfold it with a firm hand movement. To unfold the cane, hold the handle and with one firm motion lower the cane downward to fully extend it. Before use, check the connections between the individual sections of the cane.

WARNINGS:

The cane is not intended for weight-bearing support.
The product is intended for use by one person only. Operate with one hand.
The product is intended for indoor and outdoor use.

ATTENTION:

Do not hang or attach any objects to the cane. Failure to comply may result in bodily injury (e.g., during a fall) and/or damage to the product.

MAINTENANCE:

Clean the cane structure with a damp cloth using mild detergents.
Regularly remove dirt from the moving parts of the cane (height adjustment mechanism).
Do not use abrasive agents.
Do not clean the cane using pressure washers or steam cleaners.
Attention: Do not use gasoline-based agents, solvents, acids, chlorine, etc. for cleaning.
To maintain the functional properties of the product, it is important to keep all components clean and store the product properly. Storing the product in humid environments (bathrooms, laundries) or outdoors may result in loss of functional and aesthetic properties.

ATTENTION:

In the event of a "serious incident" related to the device that directly or indirectly led, could have led, or may lead to any of the following:
a) death of a patient, user, or another person, or
b) temporary or permanent serious deterioration of a patient's, user's, or another person's health, or
c) a serious public health threat
such a "serious incident" must be reported to the manufacturer and to the competent authority of the Member State where the user or patient resides.

ATTENTION:

If pain, allergic reactions, or any other disturbing or unclear symptoms occur in connection with the use of the medical device, consult a healthcare professional.

DISPOSAL AFTER WITHDRAWAL FROM USE

After the product is withdrawn from use, the medical device may be disposed of as regular municipal waste, except for electrical products – these must be disposed of in accordance with regulations for electrical and electronic equipment waste.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR MEDICAL Sp. z o.o. ,03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.
10. The warranty does not cover parts and components subject to natural wear and tear, such as wheels, tires, rubber pads, fuses, etc., the wear of which results from normal use of the product. The warranty period for batteries is 6 months.

ANTAR MEDICAL Sp. z o.o.
03-068 Warsaw, ul. Zawisłańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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