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L: 39 -45,5 cm

XL: 45- 51,5 cm

**AT04103,AT04104,AT04105 USER
MANUAL
KNEE JOINT ORTHOSIS
WITH ADJUSTABLE RANGE OF MOTION**

GENERAL WARNING

Do not use the product before and understanding these instructions. If you do not understand the warnings, cautions, or instructions, your health care professional or dealer to avoid damage to the product. Knee joint orthosis with timer made made of breathable and stretchable fabric strengthens the joint, provides compression, stabilization and gently warms the knee joint. Velcro fastening guarantees s effective immobilization. Side reinforcements with hinges ensure complete stabilization of the knee. The product has a wide range of adjustment. The orthosis can be quickly and easily adjusted to the patient's body shape.

CAUTION: Inspect the product for damage caused in transit. If you notice such damage, DO NOT use the product. For more information, contact the manufacturer.

NOTE: The manufacturer is not responsible for damage caused by negligent maintenance or cleaning, inadequate , or resulting from failure to follow the recommendations in this manual.

NOTE: Any serious incident, related to the device, should be reported to the manufacturer and the competent authority of the Member State where the user or patient resides.

CAUTION: It is forbidden to use the product other than for its intended purpose. **CAUTION:**

In the event of a product-related "serious " that directly or indirectly led to, may have led to, or may lead to any of the following:

- a) death of a patient, user or other person or
- b) temporary or permanent deterioration in the health of a patient, user or other person or
- c) a serious threat to public health

the above "serious incident" should be reported to the manufacturer and to the competent authority of the Member State in which the user or patient resides. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products. **NOTES:**

If you experience pain, allergic reactions, or other disturbing symptoms that are unclear to related to the use of the medical device, consult a healthcare professional.

CAUTION: Flammable product.

STORAGE Store the product in a dry cool room, away from moisture and direct sunlight.

TARGET GROUP OF PATIENTS Persons struggling with diseases, dysfunctions or injuries, for the treatment/rehabilitation or compensation of which this is intended (see the section on the purpose of the product of this manual). The product may be purchased by the user independently or on the recommendation of a doctor, therapist or other specialist. Whether purchasing the product on its own or as indicated by a doctor/therapist/other specialist, it is important to consider the available sizes/necessary functions and variants of the device, indications and contraindications for use, as well as the information provided by the manufacturer.

UTILIZATION

Please dispose of municipal waste in accordance with applicable general disposal conditions.

AT04105/AT04103

Indications:

- degenerative disease
 - Instability following injury and ligament damage
 - A protective measure in case of knee weakness,
 - slight sprains of the knee joint,
 - post-traumatic changes,
 - During the rehabilitation period after surgical procedures
- Contraindications: do not use the product if warming fluids or ointments have been applied, do not use in case of skin injuries.

The orthosis is equipped with two splints with adjustable flexion and extension angles, which ensure the versatility of the stabilizer for the right and left limb. Four inflexible straps help the orthosis fit the knee joint as well as possible.

Notes:

- Before putting on the orthosis, read the contents of this manual.
- Before each use, make sure that there no abrasion of the stitching and individual components on the orthosis.
- Use of the product contrary to the rules contained in the instructions can result in serious complications.
- The orthosis should be used as recommended by a doctor or physiotherapist.
- The product should not be used for open wounds, abrasions of the epidermis, or after prior application of a warming gel or ointment.
- Wearing the orthosis at night is possible only the attending physician so recommends.
- Unless otherwise advised by a doctor, the orthosis should be worn for no more than 4 hours a day.

Clock

The orthosis AT04105 has a timer with the function of adjusting the angle of flexion, extension and the possibility of complete blockade or full range of motion with metal pins. The first adjustment should take place in the presence of a specialist.

- Flexion angle adjustment in the range: 0 - 90 °.
- Upright adjustment in the range: 0 - 90 °.

Instructions for cleaning (maintenance) of the product:

- The orthosis can only be washed by hand in lukewarm water (at temperatures up to 40° C) with a mild detergent (such as gray soap).
- Remove the nibs before washing.
- Do not iron.
- Do not dry wash.
- Do not use bleach.
- Protect from moisture.
- Do not expose high temperatures.
- Do not tumble dry
- Dry in unfolded form, squeeze gently before drying.

Usage

The orthosis is worn from back, unzips completely for easy application, and then fastens each Velcro to adjust the orthosis. Length of the orthosis: 35 cm

Sizes AT04105 : (To find the right size, measure the circumference of the knee (cm). It is recommended to adjust the size by trying on the orthosis)

S: 32-35 cm

M: 35-38 cm

L: 38-42 cm

XL: 42-46 cm

XXL: 46-50 cm

XXXL: 50-54 cm

Sizes AT04103 : (To find the right size, measure the circumference of the knee (cm). It is recommended to adjust the size by trying on the orthosis)

S: 32 - 35 cm

M: 35 - 38 cm

L: 38 - 42 cm

XL: 42 - 46 cm

XXL: 46 - 50 cm

XXXL: 50 - 54 cm

Scope of delivery: orthosis, instruction manual

AT04104 Knee joint orthosis

Knee joint orthosis with clock and silicone patella reinforcement made of airy and stretchable fabric strengthens the joint, provides compression, stabilization and gently warms the knee joint. Velcro fastening guarantees effective immobilization. Side reinforcements with hinges ensure complete stabilization of the knee. The product has a wide range of adjustment. The orthosis can be quickly and easily adjusted to the patient's body shape.

Indications:

- knee injuries,
- Severe instability of the knee,
- ligament, patella conditions,
- Rehabilitation after various disorders of the knee joint,
- A protective measure in case of knee weakness,
- slight sprains of the knee joint,
- post-traumatic changes,
- During the rehabilitation period after surgical procedures.

Contraindications:

- Do not use the product if warming liquids or ointments have been applied, do not use in case of skin damage. Features:
- The orthosis ensures that the joint functions within a controlled, optimal, anatomical range of mobility,
- Side reinforcements provide additional stabilization, relieving pressure on the joint,
- Silicone reinforcement provides stabilization of the patella and prevents its pathological deformation. Notes:
- Before putting on the orthosis, read the contents of this manual.
- Before each use, sure that there is no abrasion of the stitching and individual components on the orthosis.
- Use of the product contrary to the rules contained in the instructions can result in serious complications.
- The orthosis should be used as recommended by a doctor or physiotherapist.
- The product should not be used for open wounds, abrasions of the epidermis, after prior application of a warming gel or ointment.
- Wearing the orthosis at night is possible only the attending physician so recommends.
- The duration of wearing the cortex should be consulted with the attending physician.

Clock

The orthosis has a timer with the function of adjusting the angle of flexion, extension and the possibility of complete blockade or full range of motion with metal pins. The first adjustment should take place in the presence of a specialist.

- Bend angle adjustment in the range: 0° - 90°.
- Upright adjustment in the range: 0° - 90°.

Usage

It is necessary to unfasten all the Velcro straps. The orthosis is then put on by pulling it onto the leg and then fastening the Velcro straps. Length of the orthosis: 37 cm

Scope of delivery: orthosis, instruction manual

Sizes: (To choose the right size, measure the circumference of the knee joint at the center in the upright position. It is recommended to select the size by trying on the orthosis)

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiślańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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