INSTRUKCJA OBSŁUGI
ORTEZA STAWU SKOKOWEGO I STOPY AT53076,AT53079
INSTRUCTION MANUAL
ANKLE AND FOOT ORTHOSIS AT53076,AT53079
GEBRAUCHSANWEISUNG
KNÖCHEL- UND FUSSORTHESE AT53076,AT53079
NÁVOD K POUŽITÍ
ORTÉZA KOTNÍKU A CHODIDLA AT53076,NA53079
NÁVOD NA POUŽITIE
ORTÉZA ČLENKA A CHODIDLA AT53076, AT53079
MANUEL D'INSTRUCTIONS
ORTHÈSE DE CHEVILLE ET DE PIED AT53076,AT53079
HANDLEIDING
ENKEL- EN VOETORTHESE AT53076,AT53079
MANUAL DE INSTRUCCIONES
ORTESIS DE TOBILLO Y PIE A53076,A53079
MANUALE DI ISTRUZIONI
ORTESI PER CAVIGLIA E PIEDE A53076,A53079
BRUKSANVISNING
FOTLEDS- OCH FOTORTOS VID53076,VID53079

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En

Thank you for purchasing our product. Please read the instructions for use carefully. Do not use the product until you have read and fully understood these instructions. If you do not understand the warnings, cautions, and recommendations, please contact your healthcare professional or dealer to avoid damage to the product or personal injury.

NOTE: Inspect all parts of the product for damage caused during transit. If any damage is noticed, DO NOT use the product. Instead, contact the manufacturer for more information.

AT53076

DESCRIPTION AND FEATURES

The orthosis is made of lightweight neoprene material, the product provides support for a sagging foot.

The orthosis can be easily adapted to individual needs. It provides a visible improvement in the patient's gait by providing support when the foot is raised. The orthosis can also be worn with shoes due to the use of an additional plastic part instead of the lower part of the orthosis.

INDICATIONS/APPLICATION

- Rehabilitation after surgery
- Falling foot syndrome
- Poor dorsal flexion

SIZE: UNIVERSAL

AT53079

DESCRIPTION AND FEATURES

The orthosis helps to keep the foot in a dorsiflexed position when slinging, reducing the risk of tripping and falling. It is particularly useful for people with unilateral foot drop or nerve damage. The orthosis can be used comfortably with footwear. The orthosis helps to reduce pressure around the ankle joint and also to control lateral movement.

INDICATIONS/APPLICATION

- Loss of control of the dorsal flexors of the ankle joint (heel strike/flatfoot).
- Disorders that may be associated with loss of ankle joint function e.g. CVA, herniated disc, postpolio syndrome, muscular dystrophy, spinal cord injury, polyneuropathy.

SIZE (shoe size):

S: 34-37 cm

M: 38-41 cm

L: 42-45 cm

Please note: this product is available for left and right feet.

CONTRAINDICATIONS

- Fresh or unstabilized fractures may require more restrictive forms of immobilization.
- Severe swelling or blood clots: Orthosis use may increase pressure on tissues, interfering with circulation and potentially worsening the condition.
- Open wounds or infections: The orthosis may obstruct wound access, restrict airflow, and increase infection risk.
- Sensory disorders (e.g., neuropathy): Using the orthosis without proper control may cause additional skin damage.
- Advanced arterial disease: Using the orthosis could restrict blood flow to the foot, potentially leading to serious complications like necrosis.

NOTES:

In the event of a product-related 'serious incident' that directly or indirectly led to, could have led to or is likely to lead to any of the following:

- (a) the death of a patient, user or other person; or
- (b) temporary or permanent deterioration of the health of a patient, user or other person, or
- (c) a serious threat to public health

this 'serious incident' should be reported to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. For Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

NOTES:

If you experience pain, allergic reactions or other distressing, unclear symptoms related to the use of the medical device, consult a healthcare professional.

ATTENTION: it is forbidden to use the product other than for its intended purpose

ATTENTION: The manufacturer is not responsible for damage caused by negligent cleaning and maintenance, or resulting from failure to follow the instructions in this manual.

TARGET PATIENT GROUP

Persons suffering from diseases, dysfunctions or injuries for the treatment/rehabilitation or compensation of which this device is intended (see section on intended use of the device in these instructions). The device can be purchased by the user on their own or on the recommendation of a doctor, therapist or other specialist. Whether you purchase the device yourself or on the recommendation of a doctor/therapist/other specialist, you must take into account the available sizes/necessary functions/variants of the device, the indications and contraindications for use, as well as the information provided by the manufacturer.

USE AT53076

1. Place the longer part of the orthosis on the calf and adjust to secure this part of the body firmly.
2. Using the Velcro straps, fasten the orthosis so that it does not slip off the leg or become too tight. If pain is felt, readjust the tightness using the Velcro straps.
3. Put the shorter part of the orthosis on the foot and adjust to secure this part of the body stably.
4. Use the Velcro straps to fasten the brace so that it does not slip off the leg or become too tight. If pain is felt, readjust the tightness using the Velcro straps.

USE AT53079

1. Attach the plastic part of the orthosis to the leg.
2. Adjust the leg brace with the Velcro-fastened straps so that it is not too loose and does not squeeze the leg too tightly.
3. If pain is felt, readjust the clamping force using the Velcro straps.

ADDITIONAL USE IN SHOES AT53079

1. To use the orthosis with shoes, detach the foot-attached part from the orthosis.
2. Then attach an additional flat piece to the orthosis.
3. Slide the flat accessory between the tongue of the shoe and the laces.

CLEANING AND MAINTENANCE

The product is recommended to be washed by hand with a mild detergent, at a temperature of up to 40°C. Do not wash in a washing machine. Do not use strong detergents, bleach or rough objects for cleaning. Dry laid flat in the shade, do not expose to direct sunlight. Do not tumble dry. Do not dry clean. Do not iron.

	Hand wash at up to 40 degrees C
	Do not bleach
	Do not iron
	Do not tumble dry
	Do not dry clean

STORAGE

Store the product in a dry cool room, away from moisture and direct sunlight. Not suitable for children under three years of age. Store away from sources of fire.

HOW TO DISPOSE OF THE PRODUCT

Once the device has been taken out of service, the medical device can be disposed of as normal municipal waste.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiślańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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