USER MANUAL AT02020 Back Support
AT02020 Rückenstütze GEBRAUCHSANWEISUNG
MANUEL DE L'UTILISATEUR
SOUTIEN DU DOS AT02020
AT02020 BACK SUPPORT ISTRUZIONI VAN GEBRUIK
MANUALE D'USO
SUPPORTO SCHIENA AT02020
MANUAL DE USUARIO
SOPORTE DE ESPALDA AT02020
ANVÄNDARMANUAL
BAKSTÖD AT02020

Thank you for purchasing the rear mount we offer. Before you start using the product, familiarize yourself with the following instructions for use.

The special ANTAR back support (AT02020) has been developed for people who, for various reasons, need to sit or semi-lie down for a long time and thus reduce care. (e.B. Food and drink) and discomfort is relieved. Wearing a back support can also be helpful for „healthy“ people to maintain a comfortable position, e.B. to read.

Antar rear bracket allows you to adjust the angle of inclination, is covered with a durable fabric and a special construction ensures a feeling of comfort. The back stand allows you to eat comfortably in bed, watch TV, read books or the press. It can be used by adults and older children.

The backrest can be folded completely flat, so it can be easily stored under the bed, for example.

Protect the rear kickstand from open lights and direct sunlight.

Protect the back bracket from sharp objects and surfaces.

The rear support is placed on the bed, at the headboard.

Set an angle:

Slightly lift the back bracket into the transverse tube and push slightly forward to the desired angle. Place the back foot of the stand on the desired lock. The rear bracket should now be stable.

Folds of the rear bracket:

Lift the back bracket slightly on the transverse tube and therefore release the lock. Fold the back foot over the backrest.

Cleaning instructions

Wipe the pillow with a towel. Do not use detergents or other solvents. Do not mechanically clean or dry in the dryer.

Care instructions

The product does not require maintenance.

Cut:

Backrest: 50 x 58 cm.

Base: 43 x 54 cm

WARNING:

In case a serious incident that has occurred in relation to the medical device, that directly or indirectly led, might have led or might lead to any of the following:

- a) the death of a patient, user or other person,
- b) the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- c) a serious public health threat;

the above mentioned serious incident should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established;

WARNING:

In case any pain, allergic reaction or other worrying symptoms, unclear for the user which have occurred in relation to the medical device, the user should consult a healthcare professional.

METHOD OF DISPOSAL OF THE PRODUCT AFTER REMOVING FROM SERVICE

After the device is withdrawn from use, the medical device can be disposed of as normal household waste, except for electrical products - proceed

in the manner appropriate for the disposal of electrical and electronic equipment.

ATTENTION: it is forbidden to use the product other than for the intended purpose!

THIS SIGN INDICATES THE MAXIMUM WEIGHT OF THE USER

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiańska 43, hereby grants a 12-month warranty on the rollator, starting from the date of purchase by the Buyer. Products financed by the National Health Fund are covered with the extended warranty period, which covers half of the useful life specified in the Regulation on Guaranteed Services, provided that the useful life is specified in the Regulation in months or years.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, send it immediately to ANTAR (as specified in point 1) or contact the point of sale.
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum permissible weight of the user as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts subject to normal wear and tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED. IN ORDER TO EXTEND WARRANTY PERIOD FOR PRODUCTS FINANCED BY THE NATIONAL HEALTH FUND, ATTACH A COPY OF THE ORDER FOR THE SUPPLY OF MEDICAL DEVICES BEING ORTHOPAEDIC ITEMS AND AIDS.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warszawa, ul. Zawiańska 43
Polska
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net e-mail: antar@antar.net

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v1-1 1.04.2022

Date of issue of the manual: April 1, 2022

v1-1 04/01/2022