INSTRUKCJA UŻYTKOWANIA
WKŁADKA REGULUJĄCA WYSOKOŚĆ PIĘTY DO ORTEZY STAWU SKOKOWEGO
INSTRUCTIONS FOR USE
HEEL HEIGHT ADJUSTMENT INSERT FOR ANKLE ORTHOSIS
GEBRAUCHSANWEISUNG
FERSENHÖHEN-AUSGLEICHSEINLAGE FÜR SPRUNGGELENKORTHESE
NÁVOD K POUŽITÍ
VLOŽKA PRO NASTAVENÍ VÝŠKY PATY DO ORTÉZY HLEZENNÍHO KLOUBU
NÁVOD NA POUŽITIE
VLOŽKA NA NASTAVENIE VÝŠKY PÄTY DO ORTÉZY ČLENKOVÉHO KÍBU
NOTICE D'UTILISATION
INSERT DE RÉGLAGE DE LA HAUTEUR DU TALON POUR ORTHÈSE DE CHEVILLE
GEBRUIKSAANWIJZING
INLEGZOOI VOOR HIELHOOGTEVERSTELLING VOOR ENKELORTHESE
INSTRUCCIONES DE USO
PLANTILLA REGULADORA DE ALTURA DEL TALÓN PARA ORTESIS DE TOBILLO
ISTRUZIONI PER L'USO
INSERTO PER LA REGOLAZIONE DELL'ALTEZZA DEL TALLONE PER ORTESI DELLA CAVIGLIA
BRUKSANVISNING
INLÄGG FÖR JUSTERING AV HÄLENS HÖJD FÖR FOTLEDSORTOS

MD

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En

Thank you for purchasing our product. Please read the instructions for use carefully. Do not use the product before reading and understanding this manual. If you do not understand the warnings, precautions, or recommendations, contact a healthcare professional or the seller to avoid product damage or personal injury. **WARNING:** All parts must be checked for damage caused during transport. If any damage is noticed, DO NOT use the product. For more information, contact the seller.

TARGET PATIENT GROUP

Persons suffering from diseases, dysfunctions, or injuries for the treatment/rehabilitation or compensation of which this product is intended (see the section on intended use in this manual). The product may be purchased independently by the user or based on the recommendation of a doctor, therapist, or other specialist. Whether purchased independently or recommended by a healthcare professional, the available sizes/required functions and product variants, indications and contraindications for use, as well as the information provided by the manufacturer, must be taken into account.

INTENDED USE:

The product is intended for temporary heel elevation and unloading of the Achilles tendon and the rear part of the foot.

INDICATIONS: It may be used in the following cases:

- after Achilles tendon injuries,
- after surgical procedures involving the Achilles tendon and ankle joint,
- during conservative treatment of Achilles tendon strain or inflammation,
- in overload conditions of the rear foot,
- to reduce tension in the Achilles tendon by controlled positioning of the foot in plantar flexion,
- as part of offloading therapy during the recovery period.

CONTRAINDICATIONS:

Physical or mental limitations (e.g. visual impairment) that prevent safe handling of the product. Paresthesia. Advanced varicose veins. Skin diseases or injuries at the site of application. Impaired lymphatic drainage.

LIST OF COMPONENTS: Heel height adjustment insert, instruction manual

PRODUCT FEATURES

Designed for use in a walker-type orthosis. Used for adjustable heel elevation. Stackable inserts (modular design). Fastened using a self-adhesive element. Stable adhesion to the inside of the orthosis after pressing firmly.

MODEL 1: Size: universal, left/right Adjustment range: up to 5.4 cm Material: EVA, PP

MODEL 2: Size: Universal Adjustment range: up to 5 cm Material: EVA

USE – MODEL 1

Before attaching the heel pad, clean the inside of the orthosis in that area with alcohol (e.g. denatured alcohol) and allow it to dry. Place the heel pad inside the orthosis. Remove the protective film from the dot-shaped adhesive tape located under the pad. Press the pad firmly against the inside of the orthosis. Heel pads can be stacked to achieve the desired heel elevation.

3 pads = 22° elevation

2 pads = 16° elevation

1 pad = 10° elevation

USE – MODEL 2

Before use, wash the foot and the Achilles tendon area with warm water. Dry the skin thoroughly to ensure that sweat or dirt does not affect adhesion. Position the center of the pad at the point of greatest pain or friction around the Achilles tendon. Gently press the edges to ensure full adhesion to the skin. After application, move the foot slightly to make sure the pad does not shift or curl, ensuring stable support. It is recommended not to wear the pad for more than 6 hours at a time. If redness or itching occurs, discontinue use immediately. The set includes round double-sided adhesive tapes that may facilitate stacking one pad on top of another to adjust heel height.

WARNING:

In the event of a “serious incident” related to the device that directly or indirectly led, might have led, or may lead to any of the following:

- a) death of a patient, user, or other person, or
- b) temporary or permanent serious deterioration of a patient’s, user’s, or other person’s health, or
- c) a serious public health threat,

such a serious incident must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient resides. **WARNING:** If pain, allergic reactions, or other disturbing or unclear symptoms related to the use of the medical device occur, consult a healthcare professional. **WARNING:** The product must not be used for purposes other than its intended use.

CLEANING: Hand wash only. Do not dry clean. Do not bleach. Do not tumble dry. Do not iron.

STORAGE: Store the product at a temperature between -5°C and 40°C, at 35%–75% humidity, in a dry, clean, well-ventilated place protected from direct sunlight. **DISPOSAL:** After the end of its service life, dispose of the product as regular household waste. Follow local waste segregation regulations. Protect the environment.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR MEDICAL Sp. z o.o. ,03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.
10. The warranty does not cover parts and components subject to natural wear and tear, such as wheels, tires, rubber pads, fuses, etc., the wear of which results from normal use of the product. The warranty period for batteries is 6 months.

ANTAR MEDICAL Sp. z o.o.
03-068 Warsaw, ul. Zawisłańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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